

EURODENT

A woman with brown hair tied back, wearing clear safety glasses and a white turtleneck lab coat, stands in a clinical setting. She is looking towards the camera with a slight smile. To her right is a large, modern dental unit with multiple grey handpieces and control panels. The background is a clean, white wall.

September 2021

Treat yourself!

Seize the day:

Discover the unmissable offers on the full range of Eurodent dental units

NEW

ISO/OS, the oral surgery studio

NEW

ISOFLEX, the dental unit with BESPOKE ergonomics

ISO/OS: YOUR INTEGRATED SURGICAL STUDIO

NEW




ISOCART:

Equipped with 5 instrument-holders, it can be set up according to the professional's requirements. Enables distilled water to be used for instrument sprays.



EDGE:

Surgical chair operated by 4 linear actuators via a multifunction wireless foot pedal. With Trendelenburg system and side rails.



Latest generation LED lamp that guarantees cool, long-lasting and low energy light. Provides ideal lighting whilst guaranteeing excellent working conditions and comfort



Selective, precise aspiration of blood, tissue and bone fragments. Equipped with dual jars to last longer and reduce downtime.



The unit allows instruments to be immediately available in an organised, tidy and efficient way.



An ambitious project that Eurodent has designed for professionals looking for the best instruments and equipment in the oral surgery field.

ISO FLEX: MADE-TO-MEASURE ERGONOMICS

NEW



Eurodent knows the importance of ergonomics. Now it has developed it to the highest level with ISO FLEX: this new work concept allows professionals to find the optimum position around the patient, thanks to a system that adapts to different working disciplines and habits. The freedom offered by ISO FLEX encourages maximum flexibility and fluidity at work, so the best results are achieved every time.

STARTING FROM

Shipping, assembly and VAT not included.
Promotion valid until 31/12/2021

EURODENT®

Treat yourself!

ISOPLUS DENTAL UNIT



Photos are indicative, description valid for promotion purposes.

ESSENTIAL PROMO

DENTAL CONSOLE

- Independent power and spray adjustment for each instrument.
- Safety system for automatic instrument selection.
- Instrument Anti Retraction Valve device.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes,
 - Checking of all instrument functions.
- Silicone instrument rest and removable and autoclavable handles.
- Single tray holder hanger bar on articulated arm.

INSTRUMENTATION

- Module with air - water - spray syringe.
- Module for fibre optic turbine.
- Module with EU 033 induction electric micromotor.
- P-E LED scaler.

ASSISTANT CONSOLE

- Double articulated arm.
- Two surgical suction cannulas.
- Housing for supplementary instrument.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes.
- Removable and autoclavable handle.

WATER UNIT

- Removable and autoclavable ceramic basin rotating through 160°.
- Removable and autoclavable glass filling and basin washing fountains.
- Suction preset for wet or liquid ring.

CHAIR AG04

- Electromechanical chair complete with white left armrest.
- 4 guided movements (seat up/down, backrest up/down)
- 4 programmable positions.
- Trendelenburg system synchronised with backrest.
- Double articulated headrest.
- Anti-crushing safety devices for base and back.
- Stability Kit.

MULTIFUNCTION FOOT PEDAL

- Activation and adjustment of instrument speed or power.
- Spray and chip air activation.
- Activation of chair movements and programmes.
- Activation of automatic rinse and "Last position" programmes.

EU MAIA LED SCIALYTIC LAMP

- 35000 lux light intensity, 5000° K colour temperature.
- Removable and autoclavable handles.
- Rotation on 2 axes.

E N J O Y
50%
O F F

EU MAIA LED LAMP

LED technology and a perfectly balanced rotation system enable optimum light beam positioning. Accurate definition of light spot eliminates dazzling, protecting your patients' eyes.

- Light intensity from 3000 to 35000 lux.
- Colour temperature 5000° K.
- Light spot dimensions 85x170 mm.
- Rotation on 2 axes.
- Adjustable light intensity.
- With parabolic reflector.
- Removable and autoclavable handles.



COMMERCIAL
VALUE

FREE WITH
THIS OFFER

P-E LED SCALER

Piezoelectric ultrasound fibre optic LED scaler with scaler functions, perio and endo compatible with EMS inserts.



The vibration generated is absolutely linear and longitudinal with automatic frequency search depending on the insert used.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU 033 INDUCTION MOTOR

EU 033 induction micromotor for classic conservative dentistry, preparations and prophylaxis. LED lighting.



Rotation speed between 1.000 and 40.000 rpm. Torque : 3.5 Nw/cm.

MADE IN SWISS
3-YEAR WARRANTY

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

DENTAL UNIT VALUE

-50% DISCOUNT

FINAL PRICE

JUST PER MONTH

* 72 indexed monthly instalments. All Risks insurance included.
Shipping, assembly and VAT not included.
Instruction booklets, charges and conditions at provider.
Promotion valid until 31/12/2021

EURODENT®

Treat yourself!

ISOPLUS DENTAL UNIT



Photos are indicative, description valid for promotion purposes.

SURGICAL PADDING

Thanks to the characteristics of the material used, it increases the volume of the chair's surface, adapting to the patient's body shape and distributing pressure evenly.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU-H PHOTOPOLYMERIZATION LAMP

EU-H photopolymerization lamp with two programmes (orthodontic and standard) with power adjustable between 1000 and 1800 mW/cm² and wavelength between 420 nm and 480 nm

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU MAIA LED LAMP

LED technology and a perfectly balanced rotation system enable optimum light beam positioning. Accurate definition of light spot eliminates dazzling, protecting your patients' eyes.

- Light intensity from 3000 to 35000 lux.
- Colour temperature 5000° K.
- Light spot dimensions 85x170 mm.
- Rotation on 2 axes.
- Adjustable light intensity.
- With parabolic reflector.
- Removable and autoclavable handles.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

P-E LED SCALER

Piezoelectric ultrasound fibre optic LED scaler with scaler functions, perio and endo compatible with EMS inserts.

The vibration generated is absolutely linear and longitudinal with automatic frequency search depending on the insert used.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU 033 INDUCTION MOTOR

EU 033 induction micromotor for classic conservative dentistry, preparations and prophylaxis. LED lighting.

Rotation speed between 1.000 and 40.000 rpm. Torque : 3.5 Nw/cm.

MADE IN SWISS
3-YEAR WARRANTY

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

SAFESPRAY

A pressurised container enables distilled water (with the possibility of added sanitisers) to be used for hanger bar instrument sprays and assistant's syringe.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

MUST PROMO

DENTAL CONSOLE

- Independent power and spray adjustment for each instrument.
- Safety system for automatic instrument selection.
- Instrument Anti Retraction Valve device.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes,
 - Checking of all instrument functions.
- Silicone instrument rest and removable and autoclavable handles.
- Single tray holder hanger bar on articulated arm.

INSTRUMENTATION

- Module with air - water - spray syringe.
- Module for fibre optic turbine.
- Module with EU 033 induction electric micromotor.
- P-E LED scaler.
- EU-H photopolymerization lamp.

ASSISTANT CONSOLE

- Double articulated arm.
- Two surgical suction cannulas.
- Housing for supplementary instrument.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes.
- Removable and autoclavable handle.

WATER UNIT

- Removable and autoclavable ceramic basin rotating through 160°.
- Removable and autoclavable glass filling and basin washing fountains.
- Suction preset for wet or liquid ring.
- Safespray

CHAIR AG04

- Electromechanical chair complete with white left armrest.
- 4 guided movements (seat up/down, backrest up/down)
- 4 programmable positions.
- Trendelenburg system synchronised with backrest.
- Double articulated headrest.
- Anti-crushing safety devices for base and back.
- Stability Kit.
- Surgical padding.

MULTIFUNCTION FOOT PEDAL

- Activation and adjustment of instrument speed or power.
- Spray and chip air activation.
- Activation of chair movements and programmes.
- Activation of automatic rinse and "Last position" programmes.

EU MAIA LED SCIALYTIC LAMP

- 35000 lux light intensity, 5000° K colour temperature.
- Removable and autoclavable handles.
- Rotation on 2 axes.

DENTAL UNIT VALUE

-50% DISCOUNT

FINAL PRICE

JUST PER MONTH

* 72 indexed monthly instalments. All Risks insurance included.
Shipping, assembly and VAT not included.
Instruction booklets, charges and conditions at provider.
Promotion valid until 31/12/2021

EURODENT®

Treat yourself!

ISOTRON DENTAL UNIT

3
3-YEAR WARRANTY

Photos are indicative, description valid for promotion purposes.

EU-H PHOTOPOLYMERIZATION LAMP

EU-H photopolymerization lamp with two programmes (orthodontic and standard) with power adjustable between 1000 and 1800 mW/cm² and wavelength between 420 nm and 480 nm

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU MAIA LED LAMP

LED technology and a perfectly balanced rotation system enable optimum light beam positioning. Accurate definition of light spot eliminates dazzling, protecting your patients' eyes.

- Light intensity from 3000 to 35000 lux.
- Colour temperature 5000° K.
- Light spot dimensions 85x170 mm.
- Rotation on 2 axes.
- Adjustable light intensity.
- With parabolic reflector.
- Removable and autoclavable handles.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

P-E LED SCALER

Piezoelectric ultrasound fibre optic LED scaler with scaler functions, perio and endo compatible with EMS inserts.

The vibration generated is absolutely linear and longitudinal with automatic frequency search depending on the insert used.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU 093 INDUCTION MOTOR

For endodontics and surgery: rotation speed between 80 and 40.000 rpm. The first dental micromotor with dual white and blue UV LED light. Torque: 3.5 Nw/cm - removal of old restorative composites, identifying if the resin has been completely removed during drilling; - decay detection during drilling; - composite finishing/polishing and easier removal of adhesive cements; - post-operative inspection of composite filling margins.

MADE IN SWISS 3-YEAR WARRANTY

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

SAFESPRAY

A pressurised container enables distilled water (with the possibility of added sanitisers) to be used for hanger bar instrument sprays and assistant's syringe.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

ADVANCED PROMO

DENTAL CONSOLE

- Independent power and spray adjustment for each instrument.
- Safety system for automatic instrument selection.
- Instrument Anti Retraction Valve device.
- **Display for controlling all instrument parameters to your requirements.**
- **Capacitive touch screen for controlling:**
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes,
 - Checking of all instrument functions.
- Silicone instrument rest and removable and autoclavable handles.
- Single tray holder hanger bar on articulated arm.

INSTRUMENTATION

- Module with air - water - spray syringe.
- Module for fibre optic turbine.
- **Module with EU 093 induction electric micromotor with torque control and dual white and blue UV LED.**
- P-E LED scaler.
- EU-H photopolymerization lamp.

ASSISTANT CONSOLE

- Double articulated arm.
- Two surgical suction cannulas.
- Housing for supplementary instrument.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes.
- Removable and autoclavable handle.

WATER UNIT

- Removable and autoclavable ceramic basin rotating through 160°.
- Removable and autoclavable glass filling and basin washing fountains.
- Suction preset for wet or liquid ring.
- Safespray

CHAIR AGO4

- Electromechanical chair complete with white left armrest.
- 4 guided movements (seat up/down, backrest up/down)
- 4 programmable positions.
- Trendelenburg system synchronised with backrest.
- Double articulated headrest.
- Anti-crushing safety devices for base and back.
- Stability Kit.

MULTIFUNCTION WIRELESS FOOT PEDAL

- Activation and adjustment of instrument speed or power.
- Spray and chip air activation.
- Activation of chair movements and programmes.
- Activation of automatic rinse and "Last position" programmes.
- Possibility of customising controls to your requirements.

EU MAIA LED SCIALYTIC LAMP

- 35000 lux light intensity, 5000° K colour temperature.
- Removable and autoclavable handles.
- Rotation on 2 axes.

DENTAL UNIT VALUE

-50% DISCOUNT

FINAL PRICE

JUST PER MONTH

* 72 indexed monthly instalments. All Risks insurance included.
Shipping, assembly and VAT not included.
Instruction booklets, charges and conditions at provider.
Promotion valid until 31/12/2021

E N J O Y
50%
O F F

EURODENT®

Treat yourself!

ABSOLUTE DENTAL UNIT



Photos are indicative, description valid for promotion purposes.

EXCELLENT PROMO

DENTAL CONSOLE

- Independent speed and power adjustment for motors, scalers and sprays in all instruments.
- Safety system for automatic instrument selection.
- Instrument anti retraction valve device.
- **Display for controlling all instrument parameters to your requirements.**
- **Capacitive touch screen for controlling:**
 - Chair movement and activation of chair programmes.
 - Scialytic lamp on/off.
 - Timed glass filling and basin washing.
 - Activation of automatic rinse and "Last position" programmes.
 - Checking of all instrument functions.
- Silicone instrument rest and removable and autoclavable handles.
- Double tray holder hanger bar on articulated arm.

INSTRUMENTATION

- Module with air - water - spray syringe.
- Module for proportional fibre optic turbine.
- **Module with EU 093 induction electric micromotor with torque control and dual white and blue UV LED.**
- P-E LED scaler.

ASSISTANT CONSOLE

- Double articulated and pantograph arm.
- Two surgical suction cannulas.
- Membrane keypad for controlling:
- Chair movement and activation of chair programmes.
- Scialytic lamp on/off,
- Timed glass filling and basin washing.
- Activation of automatic rinse and "Last position" programmes.
- Presetting for 2 additional instruments: syringe and polymerization lamp.

WATER UNIT

- Removable and autoclavable ceramic basin rotating through 140°.
- Removable and autoclavable glass filling fountain.
- Suction preset for wet or liquid ring, complete with centralised suction valve.
- Safespray.

CHAIR AG04

- Electromechanical chair complete with white left armrest.
- 4 guided movements (seat up/down, backrest up/down)
- 4 programmable positions.
- Trendelenburg system synchronised with backrest.
- Double articulated headrest.
- Anti-crushing safety devices for base and back.

MULTIFUNCTION WIRELESS FOOT PEDAL

- Activation and adjustment of instrument speed or power.
- Spray and chip air activation.
- Activation of chair movements and programmes.
- Activation of automatic rinse and "Last position" programmes.
- Possibility of customising controls to your requirements.

EU MAIA LED SCIALYTIC LAMP

- 35000 lux light intensity, 5000° K colour temperature.
- Removable and autoclavable handles.
- Rotation on 2 axes.

E N J O Y
50%
O F F

EU MAIA LED LAMP

LED technology and a perfectly balanced rotation system enable optimum light beam positioning. Accurate definition of light spot eliminates dazzling, protecting your patients' eyes.

- Light intensity from 3000 to 35000 lux.
- Colour temperature 5000° K.
- Light spot dimensions 85x170 mm.
- Rotation on 2 axes.
- Adjustable light intensity.
- With parabolic reflector.
- Removable and autoclavable handles.



COMMERCIAL
VALUE

FREE WITH
THIS OFFER

P-E LED SCALER

Piezoelectric ultrasound fibre optic LED scaler with scaler functions, perio and endo compatible with EMS inserts.



The vibration generated is absolutely linear and longitudinal with automatic frequency search depending on the insert used.

COMMERCIAL
VALUE

FREE WITH
THIS OFFER

EU 093 INDUCTION MOTOR

For endodontics and surgery: rotation speed between 80 and 40.000 rpm. The first dental micromotor with dual white and blue UV LED light. Torque: 3.5 Nm/cm - removal of old restorative composites, identifying if the resin has been completely removed during drilling; - decay detection during drilling; - composite finishing/polishing and easier removal of adhesive cements; - post-operative inspection of composite filling margins.



COMMERCIAL
VALUE

FREE WITH
THIS OFFER

MADE IN SWISS 3-YEAR WARRANTY

SAFESPRAY

A pressurised container enables distilled water (with the possibility of added sanitisers) to be used for hanger bar instrument sprays and assistant's syringe.



COMMERCIAL
VALUE

FREE WITH
THIS OFFER

DENTAL UNIT VALUE

-50% DISCOUNT

FINAL PRICE

JUST PER MONTH

* 72 indexed monthly instalments. All Risks insurance included.
Shipping, assembly and VAT not included.
Instruction booklets, charges and conditions at provider.
Promotion valid until 31/12/2021

EURODENT[®]

Treat yourself!

ISOPLUS DENTAL UNIT



Photos are indicative, description valid for promotion purposes.

DIGITAL 4.0 PROMO

DENTAL CONSOLE

- Independent power and spray adjustment for each instrument.
- Safety system for automatic instrument selection.
- Instrument Anti Retraction Valve device.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes,
 - Checking of all instrument functions.
- Silicone instrument rest and removable and autoclavable handles.
- Single tray holder hanger bar on articulated arm.

INSTRUMENTATION

- Module with air - water - spray syringe.
- Module for fibre optic turbine.
- Module with EU 033 induction electric micromotor.
- P-E LED scaler.

ASSISTANT CONSOLE

- Double articulated arm.
- Two surgical suction cannulas.
- Housing for supplementary instrument.
- Membrane keypad for controlling:
 - Chair movement and activation of chair programmes,
 - Scialytic lamp on/off,
 - Timed glass filling and basin washing,
 - Activation of automatic rinse and "Last position" programmes.
- Removable and autoclavable handle.

WATER UNIT

- Removable and autoclavable ceramic basin rotating through 160°.
- Removable and autoclavable glass filling and basin washing fountains.
- Suction preset for wet or liquid ring.

CHAIR AG04

- Electromechanical chair complete with white left armrest.
- 4 guided movements (seat up/down, backrest up/down)
- 4 programmable positions.
- Trendelenburg system synchronised with backrest.
- Double articulated headrest.
- Anti-crushing safety devices for base and back.
- Stability Kit.

MULTIFUNCTION FOOT PEDAL

- Activation and adjustment of instrument speed or power.
- Spray and chip air activation.
- Activation of chair movements and programmes.
- Activation of automatic rinse and "Last position" programmes.

EU MAIA LED SCIALYTIC LAMP

- 35000 lux light intensity, 5000° K colour temperature.
- Removable and autoclavable handles.
- Rotation on 2 axes.

DIGITAL 4.0 BUNDLE**

- Built-in EU Medical PC
- EU Intraoral Camera

**** ISSUE OF TECHNICAL SHEET AND SWORN APPRAISAL FOR DOCUMENTATION 4.0**

DENTAL UNIT VALUE €

-50% DISCOUNT

FINAL PRICE €

JUST PER MONTH €

* 72 indexed monthly instalments. All Risks insurance included.
Shipping, assembly and VAT not included.
Instruction booklets, charges and conditions at provider.
Promotion valid until 31/12/2021

EU ALYA-CAM LAMP WITH CAMERA

- Alya operating lamp with camera and Theia Tech, LED lamp built into the arm offering adequate lighting over the entire work area.
- LED lamp with intensity adjustable between 3000 and 50000 lux and CRI 97 for real colour perception.
- Easy, user-friendly native software for Windows, Mac and Android.
- Video and photographs shared with entire team via e-mail
 - Possibility of sharing, reviewing and making accessible various reports on main types of device
 - Primary framing directly over the area of interest (oral cavity)
 - Adequate zooming for work area.



COMMERCIAL
VALUE

BUILT-IN EU MEDICAL PC

Built-in 18.5" PC with touch screen monitor, 1366x768 pixel resolution, VESA 100x100mm, Intel I3 Dual core processor, SSD 128GB hard drive, 4GB RAM, LAN 10/100/1000 MBPS, Wireless LAN, separate 12V power unit, Windows 10 Professional, **built-in UVC LED sterilisation system.** Waterproof front certification IP65 + EN606012:2015 93/42/EEC.



COMMERCIAL
VALUE

**FREE WITH
THIS OFFER**

EU INTRAORAL CAMERA

USB digital intraoral video camera, liquid lens CMOS sensor. VGA, HD and FULL HD resolution. LED intensity, RX function, continuous/manual pre-capture digital autofocus and two-button gyro mouse controlled from the handle.



COMMERCIAL
VALUE

Programmable for use with the most frequently used image software with advanced Twain driver

**FREE WITH
THIS OFFER**

**ADDITIONAL PRICE AS AN ALTERNATIVE
TO THE INTRAORAL CAMERA**

Complete with wireless mouse and mini keypad
Complete with 120 mm fixed support
Complete with adjustable block and PC wiring

BIOS HEADREST: OPTIMUM WORKING POSITION, PATIENT COMFORT



With its special mobility and locking characteristics, **the Bios headrest** (optional) **allows you to correctly and comfortably position the patient's head, enabling maximum comfort.**

A correct working position and a comfortable, relaxed patient allow you to work with greater peace of mind and confidence.

The Bios pneumatic device enables you to swiftly adapt and lock the headrest in the optimum position with minimal effort and no mechanical tighteners. Bios follows the orbital articulations of the vertebra.

Bios achieves any working set-up and position precisely and stably.

MADE TO MEASURE FOR ME

"I chose Eurodent because it is the dental unit that allows me to fully exploit my abilities and my experience. I can customise the set-up, adapting it to my way of working and move around my workspace better. This way I am in the best condition to give better dental treatments."

Dr. Prof. Tiziano Di Carlo
Tutor on second level master's in Fixed Prosthodontics
Sapienza University of Rome

Your Eurodent Specialist is:

EURODENT[®]

Treat yourself!

Via Rossini, 26
40055 Villanova di Castenaso (BO), Italy
Tel. +39 051 780465 - Fax +39 051 781625
eurodent@eurodent.it - www.eurodent.it

Sales manager:
Ms. Giovanna Dozio
Cell. +39 393 8875521
g.dozio@eurodent.it

