

EURODENT

Treat yourself!



SUPREME

SUPREME

the only multidisciplinary chair for healthcare



SUPREME DEN
FOR DENTISTRY
TREATMENTS



SUPREME AES
FOR MEDICAL AESTHETIC
PROCEDURES



SUPREME OS
FOR ORAL SURGERY
PROCEDURES



ERGONOMICS AND COMFORT
For the patient and the dentist, who can work either standing or seated according to requirements.



FUNCTIONALITY AND DESIGN
A range of chairs suited to a variety of surgical requirements. Patient can be supine or prone. Lifting capacity 300 kg.



HYGIENE
All padded sections are easily removed for cleaning.




MADE IN ITALY
Designed and made by an Italian company renowned for its design focus and quality materials.

EURODENT[®]

Treat yourself!

REGISTERED DESIGN



TECHNICAL DATA

Minimum height	mm	470
Maximum height	mm	850
Backrest tilting angle	°	0° - 90°
Width	mm	600
Length	mm	1850
Net weight approx.	Kg	140
Maximum lifting capability	Kg	300
Lifting raising speed (with 100 Kg)	mm/sec	22
Trendelenburg	°	20°
Classification		I
Class		B



Via Rossini, 26
40055 Villanova di Castenaso (BO), Italy
Tel. +39 051 780465 - Fax +39 051 781625
eurodent@eurodent.it



www.eurodent.it